

NOMBRE PACIENTE:

Día	Mes	Año

		Nombre (Codigo)	Tipo de cartridge (precargado)	Sistema de Inyección	Tamaño de insición	Asfericidad	Filtro UV	Filtro Luz Azul	Diametro Óptico	Diametro Total	Rango de Poder	Tipo de óptica	Constante A	Poder Dioptría	Constantes de Interferometría	
MEDICONTUR	Trifocal	Esferico	Liberty Bi-Flex LIBERTY (677 MY) 	NO	MEDIET 2,2 MA	> 2,2 mm	Neutral	Si	Si	6 mm	13 mm	+8 al +30 D (inc 0,5 D) +31 al +35 D (inc 1 D)	Biconvex, aspheric, diffractive-refractive, apodized	118,9 (118,828*)	____ D	SRK/T A: 118,828 * Hoffer Q pACD: 5,431 * Holladay 1 SF: 1,682 * Haigis a0; a1; a2: 0,190; 0,192; 0,173 * Holladay 2 ACD; SF: 5,45; 1,69 ** Barret Universal II LF: 1,79 **
		Torico	Liberty Bi-Flex LIBERTY TORIC (677 MTY) 	NO	MEDIET 2,2 MA	> 2,2 mm	Neutral	Si	Si	6 mm	13 mm	+8 al +35D (inc 0,5 D) Cyl: 1,0 - 1,5 - 2,0 - 2,5 - 3,0 - 3,5 - 4,0 - 4,5 +10 al +35D (inc 0,5D) Cyl 5,25 - 6,0	Biconvex, aspheric, diffractive-refractive, apodized	118,9 (118,828*)	____ D ____ Cyl	SRK/T A: 118,828 * Hoffer Q pACD: 5,431 * Holladay 1 SF: 1,682 * Haigis a0; a1; a2: 0,190; 0,192; 0,173 * Holladay 2 ACD; SF: 5,45; 1,69 ** Barret Universal II LF: 1,79 **
	VISION EXTENDIDA (EDOF)	Esferico	3i FLEX ELON (877 PEY) 	SI	POB-MA	> 2,2 mm	Neutral	Si	Si	6 mm	13 mm	+8 al +30 D (inc 0,5 D) +31 al +35 D (inc 1 D)	Biconvex, aspheric, Refraction-based central Wavefront Linking zones	118,9	____ D	SRK/T A: 118,9 Hoffer Q pACD: 5,46 Holladay 1 SF: 1,7 Haigis a0; a1; a2: 1,32; 0,4; 0,1 Holladay 2 ACD; SF: 5,49; 1,73 Barret Universal II LF: 1,83
		Torico	3i FLEX ELON TORIC (877 PETY) 	SI	POB-MA	> 2,2 mm	Neutral	Si	Si	6 mm	13 mm	+10 al +30 D (inc 0,5D) Cyl: 1,0 - 1,5 - 2,25 - 3,0	Biconvex, aspheric, Refraction-based central Wavefront Linking zones	118,9	____ D ____ Cyl	SRK/T A: 118,9 Hoffer Q pACD: 5,46 Holladay 1 SF: 1,7 Haigis a0; a1; a2: 1,32; 0,4; 0,1 Holladay 2 ACD; SF: 5,49; 1,73 Barret Universal II LF: 1,83
PHYSIOL	Monofocal	Esferico	MICROPURE MONOFOCAL OPTIC MICROPURE 	SI	Physiol 123 Diop: +10 al +30 D. Medicel Accuject Diop: -10 al +9D y +31 al +35D	>2,2 mm	Asfericidad -0,11µ (lado posterior)	Si	Si	-10 al 24,5D 6 mm 25 al 35D 5,75 mm	-10 al 24,5D 11 mm 25 al 35D 10,75 mm	-10 al +9 D (inc 1 D) +10 al +30D (inc 0,5 D) +31 al +35D (inc 1 D)	Biconvexa Asferica -0,11µ	119,4	____ D	SRK/T A: 119.4 Hoffer Q pACD: 5.85 Holladay 1 SF: 2.06 Holladay 2 ACD: 5.85 Hagis a0; a1;a2: 1,70; 0,4; 0,1 Barret LF: 2.09
		Esferico	ISOPURE ISOFOCAL OPTIC ISOPURE 	SI	Physiol 123 Diop: +10 al +30 D. Medicel Accuject Diop: +31 al +35D	>2,2 mm	Negativa variable en su optica	Si	Si	+10 al 24,5D 6 mm 25 al 35D 5,75 mm	-10 al 24,5D 11 mm 25 al 35D 10,75 mm	+10 al +30D (inc 0,5 D) +31 al +35D (inc 1 D)	Diseño de superficie polinomial	119,4	____ D	SRK/T A: 119.4 Hoffer Q pACD: 5.85 Holladay 1 SF: 2.06 Holladay 2 ACD: 5.85 Hagis a0; a1;a2: 1,70; 0,4; 0,1 Barret LF: 2.09
	Trifocal	Esferico	FINEVISION HP TRIFOCAL OPTIC FINEVISION HI (POD F GF) 	NO	Medicel Accuject 2.2	>2,2 mm	Asfericidad -0,11µ (lado posterior)	Si	Si	6 mm	11,4 mm	+10 al +30D (inc 0.5D) +31 al +35D (inc 1D)	Biconvex Aspheric Trifocal	119,4	____ D	SRK/T A: 119.4 Hoffer Q pACD: 5.85 Holladay 1 SF: 2.06 Holladay 2 ACD: 5.85 Hagis a0; a1;a2: 1,70; 0,4; 0,1 Barret LF: 2.09
		Torico	FINEVISION HP TRIFOCAL OPTIC FINEVISION HI TORIC (POD FT 49P) 	NO	Medicel Accuject 2.2	>2,2 mm	Asfericidad -0,11µ (lado posterior)	Si	Si	6 mm	11,4 mm	+10 al +30D (inc 0.5D) +31 al +35D (inc 1D) Cyl: 1,0 - 1,5 - 2,25 - 3,0 - 3,75 - 4,5 - 5,25 - 6,0D	Biconvex Aspheric Trifocal	119,4	____ D ____ Cyl	SRK/T A: 119.4 Hoffer Q pACD: 5.85 Holladay 1 SF: 2.06 Holladay 2 ACD: 5.85 Hagis a0; a1;a2: 1,70; 0,4; 0,1 Barret LF: 2.09

* Optimized constants based on more than 400 patients

** Barrett Universal II and Holladay II constants were calculated with https://www.apacrs.org/barrett_universal2/ and <http://www.hicsoap.com> online calculators.